perfect
TENS


INSTRUCTIONS FOR USE


READ CAREFULLY BEFORE USE

TensCare™


QUICKSTART GUIDE


Connect the unit with the electrode pads
 Conecte el dispositivo a los parches
 Connectez l'appareil aux électrodes
 Verbinden Sie das Gerät mit den Elektroden

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

Dear Customer,

Thank you for choosing **perfect TENS**. TensCare stands for high-quality, thoroughly tested products for the applications in the areas of gentle electrotherapy, muscle toning, continence management and pain relief during labour.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,

Your TensCare Team

Contents

1. INTRODUCTION	5
2. INTENDED USE	5
3. PERFECT TENS FEATURES	5
4. HOW “TENS” WORKS	6
5. CONTRAINDICATIONS, WARNINGS & CAUTIONS	7
6. INFORMATION ABOUT THE PROGRAMME SETTINGS	9
7. PROGRAMMES	11
7.1. PRESET PROGRAMMES	11
7.2. MANUAL PROGRAMME	12
7.3. PROGRAMME USAGE	13
8. ELECTRODE PADS	13
8.1. PAD POSITIONING	13
8.2. PAD ADVICE	13
9. CONTENT	14
10. UNIT INFORMATION	15
10.1. CONTROLS & DISPLAY	15
10.2. OPERATING INSTRUCTIONS	16
11. SETTING UP AND USING THE PERFECT TENS	18
11.1. INSTALLATION OF BATTERIES	18
11.2. CONNECTING LEAD WIRES	19
11.3. PREPARING FOR SESSION	19
11.4. TREATMENT SESSION	19
11.5. AFTER YOUR TREATMENT SESSION	20
12. CLEANING & STORAGE	20
13. EMC	20
14. DISPOSAL OF WASTE ELECTRICAL AND ELECTRONIC PRODUCTS (WEEE)	20
15. ACCESSORIES	21
16. WARRANTY	21
17. TROUBLESHOOTING	23
18. GENERAL SPECIFICATION	24

SYMBOLS USED

	TYPE BF APPLIED PART: Equipment providing a degree of protection against electric shock, with isolated applied part. Indicates that this device has conductive contact with the end user.
	This symbol on the unit means "Refer to instructions for use".
	Temperature Limitation: indicates the temperature limits to which the medical device can be safely exposed.
	Lot Number: indicates the manufacturer's batch code so that the batch or lot can be identified.
	Humidity Limitation: indicates the humidity limits to which the medical device can be safely exposed.
	Serial Number: indicates the manufacturer's serial number so that a specific medical device can be identified.
	Do not dispose in household waste.
	Catalogue Number: indicates the manufacturer's catalogue number so that the device can be identified.
	Atmospheric Pressure: indicates the atmospheric limits to which the medical device can be safely exposed.
	Manufacturer Symbol
	Date of Manufacture: indicates the date which the medical device was manufactured. This is included within the serial number found on the device (usually in the battery compartment), either as "E/Year/Number" (YY/123456) or "E/Month/Year/Number" (MM/YY/123456).
	This medical device is indicated for home use.
IP22	This medical device is not water resistant and should be protected from liquids. The first number 2: Protected against access to hazardous parts with a finger, and the jointed test finger of 12 mm ø, 80 mm length, shall have adequate clearance from hazardous parts, and protected against solid foreign objects of 12.5 mm ø and greater. The second number 2: Protected against vertically falling water drops when enclosure is tilted up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is tilted at any angle up to 15° on either side of the vertical.
	CE Mark
	Medical Device

1. INTRODUCTION

Device Description & Principles of Design

The **perfect TENS** unit is specifically designed to provide safe and effective pain relief in the comfort of your own home. It can help manage the pain from a wide range of conditions.

Chronic pain: back pain, arthritis, rheumatic pain, lumbago, hip pain, osteoarthritic pain in the knee and sciatica central pain.

Musculoskeletal pain: fibromyalgia, falls, sprains, direct blows to muscles and auto accidents.

Acute pain: Dental and facial pain, neck pain, postoperative pain, migraine and neuralgias.

TENS machines can be used alone for pain relief or be combined with other treatments.

2. INTENDED USE

Perfect TENS is a medical device designed to be used in the home healthcare environment to provide symptomatic relief and management of acute or chronic localised pain and muscle pain and is suitable for use by all who can control the device and understand the instructions.

Do not use the device for any purpose other than this intended use.

Warning: Not suitable for use in children without medical supervision.

3. PERFECT TENS FEATURES

- **Dual Channel**

Two independent channels with four electrode pads to treat two areas at the same time.

- **Comfortable Stimulation**

Gentle TENS stimulation with 60 small steps of intensity (1 mA per step).

- **8 Preset Programmes**

Including Burst, Constant and Modulation (frequency and strength) modes and a HAN Stimulation programme.

- **1 Manual Programme**

Additional programme that allows for all the parameters to be manually adjusted and saved to suit personal needs.

- **Memory**

Features 3 functions: programme retention (automatically starts in the last programme used), number of uses and time of usage.

- **Open Circuit Detection**

Automatically resets the strength to zero if the electrode pads are not correctly applied.

- **Automatic Keypad Lock**

Automatic keypad locking prevents any accidental changes in settings.

- **Detachable Belt Clip**

Enables the user to clip the unit onto a belt.

- **Large LCD Screen**

Clearly shows the operation of the unit and the programme and intensity being used.

Note: The electrode pads need to be applied directly to the skin before the intensity increase buttons will operate. If the **perfect TENS** detects the unit is not properly connected to you (for example if the electrode pads or leads come loose) then the unit switches the intensity to zero. This safety feature prevents any unpleasant changes in output.


4. HOW “TENS” WORKS

T.E.N.S. stands for Transcutaneous Electrical Nerve Stimulation. T.E.N.S. stimulates your body's own natural defences against pain, namely the release of endorphins. TENS is totally safe and has been used successfully by thousands of pain sufferers.

TENS sends a gentle stimulation through the skin which works in TWO ways:

Pain Gate

It stimulates the sensory nerves, which carry touch and temperature signals. These nerves go to the same connections in the spine as the nerves carrying pain. A strong sensory signal will block the pain signal travelling up the spine to the brain. This is known as closing the “Pain Gate” and takes effect quite quickly after the unit is switched on. When the gate is open, pain messages get through to the brain and we feel pain. When the gate is closed, these pain messages are blocked and we do not feel pain.


Evidence suggests that TENS produce pain relief in a similar way to ‘rubbing the pain better’. The pain gate can be closed by activation of mechanoreceptors through ‘rubbing the skin’.

Scientifically, the pain gate works by release of chemical in the synapse at spinal level that inhibits transmission of pain signal.

Endorphin Release

At low frequency settings, and slightly stronger outputs, TENS drives the motor nerves to produce a small repetitive muscle

contraction. This is seen by the brain as exercise, and this promotes the release of endorphins - your body's own natural pain killer. The relief builds up and normally takes about 40 minutes to reach a maximum level which can last for hours after the machine is switched off.

By using TENS, you can expect to achieve a significant reduction in pain - if not complete relief from pain.

- TENS is effective for pain from a very wide range of causes.
- TENS machines can be used to help reduce pain from problems in muscles, joints and nerves.
- It can be also used for people with musculoskeletal pain such as long-term (chronic) back pain or knee joint arthritis. They are also often used for pain relief in the early stages of labour (see **perfect mamaTENS**),

particularly whilst a pregnant woman remains at home.

- TENS may also be used to treat many types of pain, such as migraine headaches, period pain and endometriosis (see **Ova+**), cystitis, sports injuries, fibromyalgia and neuralgia, plantar fasciitis, post-operative pain, TMJ disorder, diabetic neuropathy, osteo-arthritis and sometimes non-painful conditions such as travel sickness.
- You can use low frequency (<10 Hz) programmes on acupuncture points, to achieve similar effects to acupuncture.
- With neurogenic pain (caused by inflamed nerves) such as shingles and neuralgia, TENS may start by increasing the pain. We recommend that you only use TENS for these conditions under medical supervision.
- You can safely use TENS as long as it gives you pain relief. The effect may wear off after a few hours (this is called “accommodation”). If this happens, take a break of an hour or so before trying again. If you use settings that cause muscle movement for more than 40 minutes, you may experience aching muscles a few hours later.

5. CONTRAINDICATIONS, WARNINGS & CAUTIONS

In this manual:

A **Warning** is used when failure to follow the instructions may result in serious injury or death.

A **Caution** is used when failure to follow the instructions may result in a minor or moderate injury, or damage to the device or other property.

Notes are used to provide clarification or recommendation.

CONTRAINDICATIONS:

- **Do NOT use** if you have a pacemaker (or if you have a heart rhythm problem) or with any electronic medical devices. *Using this unit with electronic medical devices may cause erroneous operation of the device. Stimulation in the direct vicinity of an implanted device may affect some models. Stimulation on the front of the neck can affect your heart rate. Very strong stimulation across the chest may cause an extra heartbeat.*
- **Do NOT use** during the first three months of pregnancy. *It is not known whether TENS may affect foetal development.*
- **Do NOT use** on the abdomen in the later stages of pregnancy. *Stop using immediately if you experience unexpected contractions.*

WARNINGS:

Do NOT use when driving, operating machinery, or similar actions needing muscular control. *Loose electrode pads, damaged leads, or sudden changes in contact may cause brief involuntary muscle movements.*

Do NOT use to mask or relieve undiagnosed pain. *This may delay diagnosis of a progressive condition.*

Do NOT use if you have, in the area being treated: active or suspected cancer or undiagnosed pain with a history of cancer. *Stimulation directly through a confirmed or suspected malignancy should be avoided as it may stimulate growth and promote spread of cancer cells.*

CAUTIONS:

Caution should be used if you have a bleeding disorder as stimulation may increase blood flow to the stimulated region.

Caution should be used if you have suspected or diagnosed epilepsy as electrical stimulation may affect seizure threshold.

Caution should be observed when using the device at the same time as being connected to monitoring equipment with body worn electrode pads. It may interfere with the signals being monitored.

Caution should be used following recent surgical procedures. Stimulation may disrupt the healing process.

Caution Simultaneous connection to high frequency surgical equipment may

result in burns and damage to the stimulator.

Caution Strong electromagnetic fields (electro-surgery/ microwave cookers/ mobile phones) may affect the correct operation of this unit. If it appears to behave unusually, move it away from these devices.

Caution Do not permit use by children unable to understand the instructions or persons with cognitive disabilities, i.e.; Alzheimer's disease or dementia.

Caution Keep away from children under 5 years of age. *Long cord - risk of strangulation in infants.*

Caution should be observed when using the **perfect TENS** at high strength settings. Prolonged use at high settings may cause muscle injury or tissue inflammation. The **perfect TENS** is capable of delivering outputs in excess of 10 mA (RMS) or 10 V averaged over any period of 1 sec.

DO NOT PLACE ELECTRODE PADS:

- On skin, which does not have normal sensation. *If the skin is numb too great a strength may be used, which could result in skin inflammation.*
- On broken skin. *The electrode pads could encourage infection.*
- On the front of the neck. *This could cause the airway to close, giving breathing problems. May cause sudden drop in blood pressure (vasovagal response).*
- Over the eyes. *May affect eyesight or cause headaches.*
- Across the front of the head. *Effect on patients who have had strokes or seizures is not known. May affect your sense of balance. The effects of*

stimulation on the brain are unknown.

ELECTRODE PADS CAUTION:

Caution: Do not ignore any allergic reaction to the electrode pads: *If a skin irritation develops, stop using TENS, as this type of electrodes may not be suitable for you. Alternative electrode pads specially made for sensitive skin are available (see E-696-SS).*

Caution: Do not use this device with leads or electrode pads other than those recommended by the manufacturer. *Performance may vary from specification. Electrodes with smaller surface area may cause tissue irritation.*

Note: Electrode pads supplied have no measurable latex content.

TO KEEP YOUR DEVICE IN GOOD WORKING ORDER, OBSERVE THE FOLLOWING ADDITIONAL CAUTIONS:

Caution: Do not immerse your device in water or place it close to excessive heat such as a fireplace or radiant heater or sources of high humidity such as a nebulizer or kettle as this may cause it to cease to operate correctly.

Caution: Keep the device away from sunlight, as long-term exposure to sunlight may affect the rubber causing it to become less elastic and crack.

Caution: Keep the device away from lint and dust, as long-term exposure to lint or dust may affect the sockets or

cause the battery connector to develop a bad contact.

Caution: Temperature & Relative Humidity of storage: -20°C to +40°C, 8% to 70% R.H. Temperature & Relative Humidity of transportation: -20°C to +40°C, 8% to 70% R.H.

Caution: Do not attempt to open or modify the TENS unit. *This may affect the safe operation of the unit and will invalidate the warranty.*

Note: There are no known side effect to TENS use and long-term TENS use is not harmful.

6. INFORMATION ABOUT THE PROGRAMME SETTINGS

Each programme has its own combination of Frequency and Pulse Width settings which allow for different sensations through the electrode pads and suppress pain in different ways.

- **Frequency (measured in Hz - pulses per second)**

PAIN GATE: A high frequency of 110 Hz is good at blocking pain signals.

ENDORPHINE RELEASE: A low frequency of 4 or 10 Hz allows for the release of endorphins, the body's natural morphine-like substances.

- **Pulse Width (measured in µs - millionths of a second)**

The **perfect TENS** unit has pulse widths of 50 to 250 µs. Generally speaking, the higher the pulse width, the more "aggressive" the stimulation

feels, and eventually, if the pulse width is set high enough, it will usually elicit a muscle contraction, which is typically not the desired result with a TENS unit. However, if the pulse width is too low, the patient may not perceive the stimulation. Pulse Rate is important because different frequency settings target different nerve groups and the setting will determine if the "Gate Theory" or "Endorphin Theory" of TENS will be used.

- **Constant and Burst Modes**

Constant mode is when the sensation is continuous versus Burst mode when the sensation, as its name implies, is one of on and off. Constant mode is often used for acute pain via Pain Gate Effect whereas Burst mode is useful in chronic pain relief. Burst gives a combination of Pain Gate and Endorphin Release, but the squeezing feeling may not be as comfortable. The sensation in Burst is more of a 'grabbing', 'clawing' type and usually more by way of muscle twitching than with the high or low frequency modes. The stimulation intensity will need to be relatively high. In Constant mode, the sensation is continuous and is more of a tingling pins and needles type.

- **Modulation Modes**

Modulation is when either the frequency or pulse width sweep across range of settings. In modulation mode, the machine delivers a less regular pattern of TENS stimulation in an attempt to reduce or minimise the accommodation effects of regular, patterned stimulation. This is potentially most useful for patients who use TENS for several hours a day, if for no other

reason than accommodation occurs at a slower rate and therefore less intensity adjustment may be required. See programmes F and G.

- **HAN Stimulation**

Invented by Professor Han, it is an alternating mode or shifting between high and low frequencies; stimulation works at a high frequency and a relatively low pulse width; and a low frequency and a relatively high pulse width. This is thought to produce a better analgesic effect and could be used to relieve both acute and chronic pain.

Pr	Output	Frequency Hz	Pulse Width μ s
A	Constant	110	50
B	Constant	4	200
C	Burst	100 in 2Hz burst	200
D	Constant	10	200
E	Constant	110	200
F	Strength Modulated	110	200
G	Frequency Modulated	2-110 in 7 steps	200
H	HAN Alternating	100/2	100/200

* F - The user increases the strength to the desired level in the first 7 seconds. The strength drops down 90% over a period of 6 seconds and then gradually starts increasing again to the original setting over a period of 6 seconds.

* G - The frequency starts at 2 Hz and after every 3 minutes' shifts to the next frequency setting. After 21 minutes, it returns back to 2 Hz and starts the cycle over again.

Note: Specifications are nominal and subject to variation from the listed values.

7. PROGRAMMES

7.1. PRESET PROGRAMMES

There is not one programme for a particular condition, and the best choice varies from one person to another, even if they have the same type of pain.

Each user needs to select both the programme and the positioning of the electrode pads that is best for them. You may need to try a few positions/programmes before finding the one that suits you.

If you are using TENS for the first time, it is recommended that you start with programme A, which gives the gentlest stimulation.

The **perfect TENS** has 8 preset pain relief programmes.

P=Pain Gate. Gives relief when being used. Wears off after a few hours.

E= Endorphin Release. Builds up over 40 mins. Lasts several hours after use.

Programme A:

P - This is the gentlest programme. It is advisable to use for your first trial of TENS while you get used to the feeling which offers tingling/pins and needles sensation. Reduced pulse width is ideal

for treating nerve rich areas (tension type headache, facial pain, neck pain, post herpetic neuralgia, migraine).

Programme B:

E - This low frequency is also known as Acupuncture-Like. Use on acupuncture points. Pulsing or tapping sensation associated with muscle contractions. It is suitable to be used for chronic pains such as back pain & arthritis, but not for injured muscles.

Programme C:

P + E – This programme combines P and E, and is thought to be the most effective for radiating pain in arms and legs and deep muscular pain (Osteoarthritic pain in the knee, Sciatica central pain).

This programme uses high frequency in low frequency bursts. Some people find the changes in stimulation uncomfortable.

Programme D:

E - This programme provides soothing/relaxing effects and is ideal for treatment of nausea, travel and morning sickness; most successful with placing the electrode pads over the wrist (acupuncture point C6). See the instructions enclosed with the electrode pads.

Programme E:

P - Stronger, more penetrating version of programme A and a stronger feel. Could be used everywhere except very sensitive areas.

Programme F:

P - This programme continuously modulates the pulse width and gives a pleasant slow relaxing massage sensation. The continuous change helps prevent accommodation occurring.

Programme G:

P - This programme constantly varies the sensation, and is ideal for treating chronic pain over long period (example 5+ hours). Neck, shoulder, Elbow pain, Rheumatic pain, Lumbago, Hip pain, Osteoarthritic pain in the knee.

Programme H:

P + E – This programme alternates between low and high frequency every 3 seconds and gives 3 seconds of tapping, followed by 3 seconds of tingling sensation. The effective pain relief is provided by a combination of **P** and **E**. Some people find the changes in stimulation uncomfortable.

7.2. MANUAL PROGRAMME

Manual Programme M:

Surprisingly little research has been done on the importance of frequency.

Professor Han did research that showed that **P** increased as the frequency increased – but only tested up to 100 Hz.


And that **E** decreases as the frequency increases.

So, **E** is less effective when the frequency is >10 Hz and **P** is less effective when the frequency is <80 Hz.

Most clinical research on **P** has been done at 80 Hz or 100 Hz.

Many professional users believe that the response to TENS varies between individuals, and that you can get the best response by choosing the frequency that suits your nervous system the best. They suggest that you start at 80 Hz or 110 Hz, then try different the frequency until you find one that gives you better pain relief.

- Press the **P** key until you reach programme **M**.
- Press the **T** key. The Hz symbol will flash next to the number at the bottom centre of the screen.

- Press **CH1+** or **CH1-** to adjust the setting. You can select a frequency from 2 to 150 Hz.
- Press the **T** key again to save your change and advance to the next setting μ s. You can select a pulse width from 50 to 250 μ s.
- Press the **T** key to select **C** or **B**.
C = Constant **B** = Burst
- Finally, you can select the Treatment Timer. Continuous, 10, 20, 30, 40 or 60 minutes.
- Press **T** once more to return to start using.

See section 6 for more information on the choice of settings.

7.3. PROGRAMME USAGE

TENS can be used for as long as it is necessary. Continuous treatment is fine, but the electrode pads should be repositioned regularly (at least every 12 hours) to allow the skin to be exposed to the air.

The Endorphin Release programmes work better when the strength is high enough to cause small muscle contractions. Best results are achieved with a session duration between 20 and 40 minutes. Longer use may cause muscle ache.

8. ELECTRODE PADS

8.1. PAD POSITIONING

The electrode pads must always be used in pairs (two electrode pads on each channel), so that the signal can flow in a circuit.

- TENS works one vertebra at a time. You need to stimulate the sensory nerves that enter the spine at the same level as the nerve carrying your pain. Since you don't know exactly where your nerves are, the easiest way is to apply the electrode pads around/near the source of the pain.
- TENS activates the nerves best if it travels along the nerve rather than across it. So place one pad further from the spine than the source of the pain, and one closer.
- The nerves wrap around the limbs and torso, so you may have to try a few positions before you get the best effect.

- If the pain is in, or close to, your spine you can place one pad either side of the spine.
- You may feel more sensation in one pad than the other. This is normal – it depends on where the electrode pads are in relation to your nerves.

For areas that are difficult to reach, why not ask a friend to help you attach the electrode pads?

Always check that the unit is OFF before attaching or removing electrode pads.

Perfect TENS resets strength to zero if the pad or lead is disconnected from your body. This is to prevent sudden changes in sensation when the pad is re-connected.

Note: See examples of positioning in the instructions enclosed with your electrode pads.

8.2. PAD ADVICE

- The electrode pads supplied are reusable but for single patient use.
- In order to obtain the best conductivity through the electrode pads always ensure that they are in good condition and tacky.
- Before use make sure your skin is clean and dry.
- Peel the electrode pads from their protective plastic shield by holding and lifting one corner of the pad and pulling. **Do not** pull on the pigtail wire of the pad.

- After use always replace the electrode pads on the plastic liner and replace them in the re-sealable plastic bag.
- If the electrode pads dry out, then it is best to buy a replacement pack. In an emergency, it may be possible to restore some of the tackiness of the pad by adding a tiny drop of water on each pad and spreading around. If too much water is added, the electrode pad will become too soft. If that happens then it is suggested in order to try and re-establish some adhesiveness, to place them sticky side up in a refrigerator for a few hours.
- In very hot weather the gel on the electrode pads may become soft. In such cases place them, still on their plastic liners and in their pouch into a fridge until they return to their normal condition.
- The electrode pads provided are latex-free.
- Replace the electrode pads when they lose their stickiness.
Poor connection may cause discomfort and skin irritation.
- Storage life of an unopened pack of electrode pads is 2 years. This may be affected by very high temperatures or very low humidity.

9. CONTENT

The pack contains:

- 1 x **Perfect TENS** pain relief unit
- 2 x Lead wires (L-CPT)
- 4 x 50x50 mm electrode pads (E-CM5050)
- 2 x AA 1.5V alkaline batteries
- 1 x Detachable belt clip
- 1 x Storage pouch
- 1 x Manual instruction

10. UNIT INFORMATION

10.1. CONTROLS & DISPLAY

10.2. OPERATING INSTRUCTIONS

ON/OFF

To turn the unit **ON**, press the **ON/OFF** button on the top of the unit and hold for 3 to 5 seconds until the display shows.

To turn the unit **OFF**, press the **ON/OFF** button on the top of the unit and hold for 3 to 5 seconds until the display stops

The unit will turn **OFF** automatically:

- When the Timer reaches zero,
- If it is left at zero strength for more than 5 minutes.

The unit will start with the programme you used last.

At first use, or after changing batteries, the display shows that the unit is automatically set in Programme **A** at zero strength.


Note: Always check unit is **OFF** before applying or removing electrode pads.

The backlight will turn off 10 seconds after the last button press.


On the front of the unit there are six buttons:

STRENGTH CONTROLS

The buttons marked **CH1+**, **CH1-**, **CH2+**, **CH2-** are the strength controls.

The left-hand buttons change strength in the left-hand lead, CH1, and the right-hand buttons change the right-hand lead, CH2.

To increase strength in steps of 1 mA, press and release the top part of the two channel buttons, marked **CH1+** and **CH2+**.

The strength levels are shown on the LCD.

The strength control buttons will not operate until the unit is properly connected to you. Perfect TENS detects a disconnection and automatically returns the strength to zero.

The unit has 60 levels of strength. If you hold down the button for 3 to 5 seconds, the strength will start scrolling.

You may feel nothing over the first few presses. Continue pressing until the sensation is strong but comfortable. Further increases during use may be necessary if your body becomes used to the sensation.

To decrease the strength, press and release the lower part of the buttons marked **CH1-** and **CH2-**.

The yellow LED on the output socket indicates that there is an active output. The display will remain on for 5 seconds after the plug is removed.

PROGRAMME CONTROL

The button marked **P** is the programme control. The **perfect TENS** has eight different, specially developed programmes lettered **A** to **H**, plus a manual programme **M**. When the unit is first switched on it automatically enters programme **A**.

Each time you press and release the **P** button, the programme changes and the letter is shown on the LCD.

Each time you change the programme, the strength level reverts back to zero. This is a safety feature to alleviate any sudden feeling of a surge, as each programme gives a different sensation.

TREATMENT TIMER

The button marked **T** can be used to set the session duration and to select manual parameters for adjustment. When you switch the unit on, it is automatically set at Continuous, (displayed as C) in all programmes but **B, C & H**, which means that the unit will continue functioning until you switch it off.

When you press the **T** button, the clock symbol will flash.

You can set session times of 10, 20, 30, 40 or 60 minutes with the **CH1+** and **CH1-** buttons.

Press **T** again to save your selection.

The LCD shows the session duration next to the clock symbol. The unit automatically counts down the minutes set and switches off when it reaches 0.

LOW BATTERY

An empty battery symbol will show when you need to change the batteries. The unit will shut down about 2 minutes after this.

KEYPAD LOCK

If you do not press any keys for 30 seconds, the keypad will lock so that only the keys marked **CH1-** or **CH2-** will work. This is to avoid accidental changes in setting.

To unlock, press either **CH1-** or **CH2-** keys.

OPEN CIRCUIT CUTOFF

If the **perfect TENS** is not correctly connected to your body, the strength in the channel which is not properly connected will automatically reset to zero. This is to prevent sudden changes if a broken connection is re-made.

MEMORY

The **perfect TENS** has a Memory with three functions:

1) Programme Retention. When you switch the unit on, it will automatically start in the programme which was being used when it was switched off.

2) Usage. Press **T** and **CH2-** together and hold down for 3 to 5 seconds. The display will show the number of times used and duration in hours used.

Press the same buttons again to return to normal controls.

3) Memory Reset. To reset memory to zero, hold down the **T** and **CH1-** buttons together for 3 to 5 seconds.

11. SETTING UP AND USING THE PERFECT TENS

11.1. INSTALLATION OF BATTERIES


- 1) Remove belt clip by sliding down.

- 2) Remove battery cover by pulling on tag.

- 3) Insert batteries.

Ensure that the batteries are inserted the right way as shown in battery compartment and that the ribbon is behind them.

- 4) Replace battery cover and belt clip.

When the batteries are running low, a low battery indicator will show on the screen and it is important to change the batteries as soon as possible.

Rechargeable batteries

The unit will work with rechargeable batteries, but the display may appear dim.

Storage

Remove batteries from your **perfect TENS** if the unit is unlikely to be used for a long period. *Some types of batteries may leak corrosive fluid.*

Battery Life

Batteries should last about 48 hours at full power.

Unused batteries have a nominal shelf life of 3 years, but will usually last longer than this.

Battery Warnings

Do NOT pierce, open, disassemble, or use in a humid and/or corrosive environment.

Do NOT expose to temperatures over 60°C(140F).

Do NOT put, store or leave near sources of heat, in direct strong sunlight, in a high temperature location, in a pressurized container or in a microwave oven.

Do NOT immerse in water or sea water or get wet.

Do NOT short-circuit.

Do NOT connect the device unless the battery cover is in place.

If battery leakage occurs and comes in contact with the skin or eyes, wash thoroughly with lots of water.

Keep batteries out of the reach of small children.

Caution NEVER attempt to recharge an alkaline battery. *Risk of explosion.*

Caution Do not mix old, new or different types of batteries as this may lead to battery leakage or low battery indication.


Disposal: Always dispose of batteries responsibly according to local government guidelines. Do not throw batteries onto a fire. *Risk of explosion.*

11.2. CONNECTING LEAD WIRES

Insert the lead wire plug into the base of the unit.

Attaching the electrode pads to the lead wire.

Push the pin ends firmly into the pigtail ends of the electrode pads.

The lead wires may be damaged by rough handling, and should be treated with care.

Lead wire colour coding.

The ends of the lead wires are coloured white or red. This coding is provided for some professional uses. For most purposes, the orientation makes no difference, and you can ignore this colour coding.

11.3. PREPARING FOR SESSION

- 1) Before use make sure your skin is clean and dry.
- 2) Place the electrode pads on the skin, see section 8.
- 3) You can use your TENS for several hours, you can attach the unit to your belt using the belt clip or hang it around your neck using a neck cord (not included). For safety, use thread or wool that will snap easily if the unit snags on something.
- 4) You can safely lean on the electrode pads. The sensation might change.

Warning: Ensure the **perfect TENS** is switched OFF before applying the electrode pads on the skin.

11.4. TREATMENT SESSION

- 1) Press and hold the **ON** button on the control unit for 3 to 5 seconds to switch the control unit on.

When switched on for the first time, the **perfect TENS** will automatically select the programme **A**. Next time, it will automatically select the programme you were using the last time it was switched off.

- 2) You can select from the 8 preset programmes or the manual programme. Details in section 7 will help you identify the best programme to suit you.
- 3) With the required programme selected, you can adjust the intensity of the stimulation with the buttons **CH1-**, **CH1+**, **CH2-** and **CH2+** until you reach

a comfortable level. The LCD display shows the strength of intensity used for each channel.

Note: If the sensation becomes uncomfortable, reduce the intensity.

11.5. AFTER YOUR TREATMENT SESSION

When the timer reaches 0, your treatment session is complete and the unit turns off.

1) Check that the control unit is OFF. If it is not, hold down the **OFF** button to switch off then remove the electrode pads from your skin by holding the pad itself and gently pulling. Be careful not to pull the wires as this may damage them.

2) Replace the electrode pads to their protective plastic shield and return them to the re-sealable plastic bag.

Note: When removing the electrode pads, DO NOT PULL ON LEAD WIRES.

12. CLEANING & STORAGE

Clean the case and lead wires at least once a week by wiping with a damp cloth and a solution of mild soap and water. Wipe dry.

- **Do not** immerse your TENS machine in water.
- **Do not** use any other cleaning solution than soap and water.

Storage life

- Storage life of an unopened pack of self-adhesive electrode pads is 2 years. This may be affected by very high temperatures or very low humidity.
- The unit has no fixed shelf life.

13. EMC

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d = 3,3$ m away from the equipment.

(Note. As indicated in Table 6 of IEC 60601-1-2:2007 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields $d = 3,3$ m at an IMMUNITY LEVEL of 3 V/m).

Note: For hospital use, full EMC advice tables are available on request.

14. DISPOSAL OF WASTE ELECTRICAL AND ELECTRONIC PRODUCTS (WEEE)

One of the provisions of the European Directive 2002/96/CE is that anything electrical or electronic should not be treated as domestic waste and simply thrown away. To remind you of this Directive all affected products are now being marked with a crossed-out wheellie bin symbol, as depicted below.

To comply with the Directive, you can return your old electro-therapy unit to us for disposal. Simply print a postage-paid PACKETPOST RETURNS label from our website www.tenscare.co.uk, attach this to an envelope or padded bag with the unit enclosed, and post it back to us. Upon receipt, we will process your old device for components recovery and recycling to help conserve the world's resources and minimise adverse effects on the environment.

15. ACCESSORIES

Expected Service Life

- The machine will often last for more than 5 years, but is warranted for 2 years. Accessories (lead wires, electrode pads, and batteries) are not covered by the warranty.
- Lead life depends greatly on use. Always handle the leads with care. We recommend to replace the lead wires regularly (about every 6 months).
- Electrode pads should last 12 to 20 applications, depending on skin condition and humidity.
- AA alkaline batteries should last about 48 hours of continuous use.

Replacement electrode pads, new batteries and lead wires are available from your supplier or distributor (see back cover for contact details), by mail order from TensCare, by telephone using a credit or debit card, or through the TensCare website.

The following replacement parts may be ordered from TensCare at www.tenscare.co.uk or +44(0) 1372 723434.

E-CM5050	Pack of 4 50x50 mm electrode pads
E-696-SS	Pack of 4 50x50 mm hypoallergenic electrode pads
L-CPT	Lead wire
B-AA	1.5V AA batteries
X- BC-PT	Replacement belt clip
X- BL-PTT	Replacement battery cover

16. WARRANTY

This warranty refers to the unit only. It does not cover, electrode pads, battery, or lead wires.

PRODUCT INFORMATION	WARRANTY
---------------------	----------

This product is warranted to be free from manufacturing defects for 2 years from date of purchase.

This warranty is void if the product is modified or altered, is subject to misuse or abuse; damaged in transit; lack of responsible care; is dropped; if incorrect battery has been fitted; if the unit has been immersed in water; if damage occurs by reason of failure to

follow the written instructions for use booklet enclosed; or if product repairs are carried out without authority from TensCare Ltd.

We will repair, or at our option replace free of charge, any parts necessary to correct material or workmanship, or replace the entire unit and return to you during the period of the warranty. Otherwise, we will quote for any repair which will be carried out on acceptance of our quotation. The benefits conferred by this warranty are in addition to all other rights and remedies in respect of the product, which the consumer has under the Consumer Protection Act 1987.

Our goods come with warranties that cannot be excluded under the UK consumer Law. You are entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality.

Before you send your unit for service

Before sending in your unit for service, please take a few minutes to do the following:

Read your manual and make sure you follow all the instructions for use.

Returning your unit for service

Should repair be needed within the warranty period, enclose the tear off section of the warranty card (see page 25) and your proof of purchase receipt. Please ensure all relevant details are completed before sending your unit in for service. Please ensure your contact

details are still current and include a brief description of the problem you are experiencing together with your purchase receipt.

For hygiene reasons, please do not include used electrode pads. Send only the unit and lead wires.

Please return the unit and warranty card (see page 25) at your cost to:

TensCare Ltd

PainAway House,
9 Blenheim Road,
Longmead Business Park,
Epsom, Surrey
KT19 9BE, UK

Should you require any further information please do not hesitate to contact us by calling our number:

+44 (0) 1372 723 434.

17. TROUBLESHOOTING

If your TENS machine is not working properly, please check the following:

Problem	Possible causes	Solution
No display	Flat batteries.	Replace batteries.
	Batteries inserted incorrectly.	Remove plastic wrap. Check + / - .
	Damaged springs in battery compartment.	Contact supplier.
Low battery display	Low batteries.	Replace batteries.
No sensation	Intensity is not strong enough.	Increase strength. Most users will feel something at a setting below 20.
Output will not increase above zero	Open circuit cut-out operating. Lead not connected to body or faulty/damaged.	Ensure machine is attached to your body correctly. Try using the second lead wire. Purchase replacement if necessary.
	(Lead wires can break at the bend where they leave the machine giving no, or intermittent, output).	
Sudden change in sensation	If you disconnect and re-connect a few minutes later, the signal will feel quite a lot stronger.	Always return strength to zero after disconnecting the leads or electrode pads.

If the above review has failed to resolve your problem, or to report unexpected operation or events, call TensCare or your local supplier or distributor (address on back cover) for advice.

Contact TensCare customer service on +44 (0) 1372 723 434. Our staff are trained to assist you with most issues you may have experienced, without the need to send your product in for service.

European Medical Device Regulation requires that any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority in your country. This can be found at: <https://ec.europa.eu/docsroom/documents/36683/attachments/1/translations/en/renditions/pdf>

18. GENERAL SPECIFICATION

Waveform	Asymmetrical Bi-Phasic
Amplitude (over 500 Ohm load)	60 mA zero to peak +ve in 60 steps. +/- 10%
Output plug	Fully shielded: touch proof
Channels	Dual channel
Batteries	2 x AA alkaline (two AA batteries)
Weight	75 g without batteries
Dimensions	120 x 60 x 20 mm
Safety Classification	Internal power source. Designed for continuous use. No special moisture protection.
Environmental specifications:	
Operating:	Temperature range: 5 to 40°C Humidity: 15 to 93% RH non-condensing Atmospheric pressure: 70 kPa to 106 kPa
Storage:	Temperature range: -25 to +70°C Humidity: Up to 93% RH non-condensing
TYPE BF APPLIED PART 	Equipment providing a degree of protection against electric shock, with isolated applied part.
	This symbol on the unit means "Refer to Instructions for use".
IP22	The unit is not water resistant, and should be protected from liquids.
	Complies with EU WEEE regulations

Note: The electrical specifications are nominal and subject to variation from the listed values due to normal production tolerances of at least 5%.

PLEASE RETAIN THIS WARRANTY CARD.

RETURN THIS PORTION ONLY WHEN YOU RETURN YOUR PRODUCT FOR REPAIR UNDER WARRANTY.

NAME: _____

ADDRESS: _____

POSTCODE: _____

DAYTIME TELEPHONE: _____

E-MAIL: _____

MODEL: _____

DATE OF PURCHASE: _____

ATTACH PROOF OF PURCHASE
DO NOT SEND IN ELECTRODE PADS

RETAILERS NAME: _____

RETAILERS ADDRESS: _____

RETAILERS POSTCODE: _____

BRIEF DESCRIPTION OF PROBLEM YOU ARE EXPERIENCING: _____

WARRANTY IS VOID UNLESS THE ABOVE INFORMATION IS COMPLETED AND CORRECT.

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TensCare aim to give you the best possible product and service. We listen to your suggestions and are constantly trying to improve our products. We also want to learn about the way our products are used, and the benefits they give. If you have anything you would like to share with us, please get in touch.

www.tenscare.co.uk

Follow us:

TensCare Ltd

@TensCareLtd

TensCare Limited

@tenscareltd

EC Declaration of Conformity

TensCare Ltd hereby declare that an examination of the production quality assurance system has been carried out following the requirements of the UK national legislation according to Annex V of the Directive 93/42/EEC on medical devices. We certify that the production quality system conforms with the relevant provisions of the aforementioned legislation, and the result entitles the organization to use the CE 2797 marking on this product.

Distributed by:

TensCare Ltd

9 Blenheim Road, Epsom,
Surrey KT19 9BE, UK
Tel: +44(0) 1372 723434
www.tenscare.co.uk

Advena Ltd

Tower Business Centre,
2nd Flr, Tower Street,
Swatar, BKR 4013, Malta