



BSU895UK

MAIN POWER SUPPLY ADAPTER 14V 2P UK

- Specific main power supply adapter for BSU suction units family
- 2-pole male connection cable complete with retaining ring
- Wide range input voltage: 100÷240VAc can be used at 50 or 60 Hz
- Allows the device to be powered and the internal battery to be recharged
- It can be always connected to the medical suction unit
- Tested to comply with the latest international regulatory standards
- Rated power: 60 W. Maximum power: 75 W
- Output voltage 14 Vcc. Delivered rated power: 5 Adc
- It can be connected the suction unit while operating on the patient
- It is an accessory component of the BSU medical devices and is CE marked. The LYD power supply unit complies with IEC 60601-1-2: 2014 and the standards of the IEC 60601 family
- The user can change the cable to make it suitable for the type of sockets available in his state of use (type C8 socket)
- The power supply is not replaceable with other similar models even if of the same brand
- The device has a warranty of 12 months from the date of sale and cannot be repaired (excluding the input connector)

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.