

ALASKA XL RISER RECLINER

- **Dual Motor**

A **HUGE**
addition to our
Riser Recliner
Range!!



The Alaska Extra Large Dual Motor Riser Recliner features an extra wide seat of 72cm (28.3") and supports users up to 200kg (31.5st), providing a seating solution for people with larger physiques to relax in comfort.



LL512

www.drivedevilbiss.co.uk

drive

DeVilbiss
HEALTHCARE



The Alaska is upholstered in a wipe-down Anti-Microbial PVC, which is Crib 5 fire safety compliant, making it ideal for use in multi-user care environments.

Features:

- » The widest riser in our range with a generous seat width of 72cm (28.3")
- » The highest maximum user weight in our range, up to 200kg (31.5st)
- » Dual motor and 3-way mechanism enables the legrest and backrest to be operated independently, providing countless recline positions
- » 5 button handset with Smart Function enabling the chair to return to standard sitting position at the touch of a button.
- » Waterfall and ergonomically contoured backrest cushions
- » Provides a huge selection of comfortable recline positions for watching TV or reading through to an excellent "Snooze" position
- » Fully padded chaise for increased comfort
- » 2 Convenient side pockets for easy storage
- » Safety key provides a motion lockout function to prevent unintended movement
- » In-built battery operated reset function (Battery Back up) allows the chair to move to any position in the event of a power cut, for one time max**
- » Quiet motors provide a smooth and discreet operation
- » Upholstered in Anti-Microbial PVC, for increased hygiene.
- » Upholstery is Crib 5 fire safety compliant, making the chair ideal for use in multi-user care environments.
- » Available in Black and Cedar brown* (Cedar brown is subject to lead times)
- » 2-year warranty on frame and electrical components

**Batteries not included and require two 9V PP3 batteries

Height 115cm (45.3")	Width 103cm (40.5")	Depth 89cm (35")	Max User Weight 200kg (31.5st)
Seat Height 48cm (19")	Seat Width 72cm (28.3")	Seat Depth 54cm (21.3")	

If you would like to view this product please contact your local stockist.



2 Handy
Storage Pockets



Extra Wide
72cm (28.3") Seat



Comfortable
TV Position



Full Snooze
Position



Comfortable
Padded Upholstery



User Friendly
Hand Control



Fully Padded
Chaise



Anti-Microbial
PVC

Available in 2 AM PVC Colours



Black
Code: CLR50BLK



Cedar Brown*
Code: CLR50CDB

*Cedar brown is subject to lead times

All the dimensions are for the armchairs in their normal chair position. Tolerance is ± 25 mm. All units should be placed with the back at least 60cm away from the wall to prevent damage when reclining the chair. Your riser recliner armchair complies with The Furniture and Furnishings (Fire) (Safety) Regulations 1988.