

## Technical Specifications

- Temperature sensor: semiconductor (0-45°C)
- Flow sensor: bi-directional digital turbine
- Flow range:  $\pm 16$  L/s
- Volume accuracy:  $\pm 3\%$  or 50 mL
- Flow accuracy:  $\pm 5\%$  or 200 mL/s
- Dynamic resistance at 12 L/s:  $<0.5$  cm H<sub>2</sub>O/L/sec
- Communication port: USB
- Power Supply: line powered from USB port
- Dimension: 50 x 142 x 26 mm
- Weight: 65 g