

Technical Specifications

Stainless steel chamber: Ø 24.5 x 32 cm

External size: 51 x 39 x 59 cm (W x H x D)

Net weight: 54 Kg

2 sterilization cycles: 121°C and 134°C

Operating Voltage: 230 V ± 10% - 50/60 Hz

Power consumption: 1500 W

Class N

Made in Italy

Warranty: 1 year

CYCLES NOT WRAPPED 121°C

Materials:

Fragile solids

Rubber solids.

Glass

Turbines if indicated by producer

Sterilization time: 35 min.

Max Load: 4 Kg (1 Kg load in any tray)

Working pressure: 1.2 ÷ 1.3 bars

CYCLES NOT WRAPPED 134°C

Materials:

Stainless metal solids

Turbines if indicated by producer

Sterilization time: 21 min.

Max Load: 4 Kg (1 Kg load in any tray)

Working pressure: 2.1 ÷ 2.2 bars

