


# Bella lumax



Bella lumax  
on/off switch

Maximum power in a lamp format examination lamp. Same power as the Bellon model with the advantage of fitting all LID fixtures.

The 50/100% intensity option allows you to switch from examination light to intervention light

● Medical device class 1

Light intensity **2300 Lumens\***

Colour of temperature **4500°K\* (neutral white)**


Colour rendering index **CRI 95/100**

Electric consumption **18 Watts**

- LED life span: 50 000 h minimum
- Power cable length: 3 meters / 4 meters
- Suggested colors: white
- Head dimensions: 16 x 13 cm
- Electrical protection: insulation class II
- Standards: IEC 60 601 -1/-2-41

- Fixings to choose from p23-26



| Lux*    | Beam 10°                                                                            | Distance | Illuminated diameter |
|---------|-------------------------------------------------------------------------------------|----------|----------------------|
| 240 000 |  | 30 cm    | 10 cm                |
| 160 000 |                                                                                     | 50 cm    | 13 cm                |
| 44 000  |                                                                                     | 100 cm   | 20 cm                |

\* Values accurate within +/- 10 %

| CODES                         |               |                    |
|-------------------------------|---------------|--------------------|
|                               | On/off Switch | Contactless switch |
| BELLA LUMAX<br>Flexible 70 cm | LED22650      | LED22650-ISC       |
| + Double intensity option     | LED22650DI    | LED22650DI-ISC     |

