

LUXIFLEX SENSOR PLUS

Luxiflex Sensor Plus is a cold light luminaire, versatile and with a compact design. Luxiflex Sensor Plus is manufactured with a contactless switch. This luminaire is designed for examination procedures that require great performance.

Flexibility

Flexible arm allows to adjust the height to assure the proper working distance.

Compact design

Luxiflex offers easy use, thanks to its compact LED head with a diameter of 8 cm and its contactless switch located at the base of the lamp. These features make this model an ideal luminaire for the most demanding examinations procedures.

Mounting options

The standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolleys stand and rail brackets.

High performance LED

Luxiflex Sensor Plus 35.000 lux at 50 centimeters.

Contactless switch

Allows to switch the lamp on and off without contact.



LUXIFLEX PLUS SENSOR

Technical Details (Cod. I0104NLSR)	
Light source	LED 10W 14°
Lamp lifetime	50.000 hours
Lux a 50 centimeters	35.000 lux
Colour temperature	4.500 °K
Colour rendering	97%
Total length	85 cm
Colour	White



Contactless switch



Trolley stand
4 kg Cod. 09655
8,8 kg Cod. 09651



Wall support
 Cod. B0101800



AH table clamp
 Cod. B010190



Rail Plus bracket
 Cod. 82026



Extension wall arm
 Cod. B0101S0